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1. クオンティフェロン TB-2Gの使用指針, 結核. 2006; 81:393-7.
2. Harada N, Nakajima Y, Higuchi K, Sekiya Y, Rothel J, Mori T. Screening for tuberculosis infection using whole-blood interferon-gamma and Mantoux testing among Japanese healthcare workers. *Infect Control Hosp Epidemiol* 2006; 27:442-8.
3. 大森正子, 星野齊之, 山内祐子, 内村和広. 職場の結核の疫学的動向—看護師の結核発病リスクの検討—結核. 2007; 82:85-93.
4. Hotta K, Ogura T, Nishii K, Kodani T, Onishi M, Shimizu Y, Kanehiro A, Kiura K, Tanimoto M, Tobe K. Whole blood interferon-gamma assay for baseline tuberculosis screening among Japanese healthcare students. *PLoS ONE* 2007; 2:803.
5. Farhat M, Greenaway C, Pai M, Menzies D. False-positive tuberculin skin tests: what is the absolute effect of BCG and non-tuberculous mycobacteria? *Int J Tuberc Lung Dis* 2006; 10:1192-204.
6. Harboe M, Oettinger T, Wiker HG, Rosenkrands I, Andersen P. Evidence for occurrence of the ESAT-6 protein in *Mycobacterium tuberculosis* and virulent *Mycobacterium bovis* and for its absence in *Mycobacterium bovis* BCG. *Infect Immun* 1996; 64:16-22.
7. Mori T, Sakatani M, Yamagishi F, Takashima T, Kawabe Y, Nagao K, Shigeto E, Harada N, Mitarai S, Okada M, Suzuki K, Inoue Y, Tsuyuguchi K, Sasaki Y, Mazurek GH, Tsuyuguchi I. Specific detection of tuberculosis infection: an interferon-gamma based assay using new antigens. *Am J Respir Crit Care Med* 2004; 170:59-64.
8. Andersen P, Munk ME, Pollock JM, Doherty TM. Specific immune-based diagnosis of tuberculosis. *Lancet* 2000; 356:1099-1104.
9. Joseph HA, Shrestha-Kuwahara R, Lowry D, Lambert LA, Panlilio AL, Raucher BG, Holcombe JM, Poujade J, Rasmussen DM, Wilce M. Factors influencing health care workers' adherence to work site tuberculosis screening and treatment policies. *Am J Infect Control*. 2004; 32:456-61.
10. Kendig EL Jr, Kirkpatrick BV, Carter WH, Hill FA, Caldwell K, Entwistle M. Underreading of the tuberculin skin test reaction. *Chest* 1998; 113:1175-7.
11. QuantiFERON®-TB Gold Package Insert, reviewed and approved by the FDA.
12. Diel R, Lodenkemper R, Meywald-Walter K, Niemann S, Nienhaus A. Predictive Value of a Whole-blood IFN- γ Assay for the Development of Active TB Disease. *Am J Respir Crit Care Med* 2008; 177:1164-70.
13. Nienhaus A, Schablon A, Le Bâcle C, Siano B, Diel R. Evaluation of the interferon-gamma release assay in healthcare workers. *Int Arch Occup Environ Health* 2008; 81:295-300.
14. Mazurek GH, Jereb J, Lobue P, Iademarco MF, Metchock B, Vernon A. Guidelines for using the QuantiFERON-TB Gold test for detecting *Mycobacterium tuberculosis* infection, United States. *MMWR Recomm Rep* 2005; 54 (RR-15):49-55.
15. Berrington de González A, Darby S. Risk of cancer from diagnostic X-rays: estimates for the UK and 14 other countries. *Lancet*. 2004; 363:345-51.

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