

# Analysis Software\* Quick Reference Guide

Run Date: Tuesday, 1 January 2008  
 Run Number: Assay 1  
 Kit Batch Number: FFS400401  
 Operator: FD A

Enter test information.

Paste Raw Data

Copy and paste data, or enter values manually.

|   | 1     | 2     | 3     | 4     | 5     | 6     | 7     | 8     |
|---|-------|-------|-------|-------|-------|-------|-------|-------|
| A | 0.065 | 0.103 | 3.012 | 1.864 | 1.999 | 0.071 | 0.256 | 0.166 |
| B | 2.103 | 2.314 | 2.997 | 1.370 | 1.352 | 0.102 | 0.400 | 0.253 |
| C | 0.053 | 0.190 | 2.460 | 0.913 | 0.937 | 0.057 | 0.061 | 1.432 |
| D | 0.052 | 0.193 | 1.379 | 0.512 | 0.507 | 0.062 | 0.062 | 3.009 |
| E | 0.053 | 0.079 | 3.005 | 0.283 | 0.231 | 0.052 | 0.195 | 3.001 |
| F | 0.065 | 0.167 | 2.999 | 0.109 | 0.154 | 0.937 | 1.023 | 0.951 |
| G | 0.074 | 0.097 | 0.222 | 0.089 | 0.102 | 0.075 | 0.079 | 0.262 |

Create and save format (default or manual method).

Change Subject ID  
 Subject ID: A0090005  
 OK Cancel

Enter subject names. Calculate results.

| Std | Conc  | Mean  | % CV | QC Result |
|-----|-------|-------|------|-----------|
| S1  | 8.00  | 1.932 | 4.9  | PASS      |
| S2  | 4.00  | 1.366 | 0.4  | PASS      |
| S3  | 2.00  | 0.925 | 1.8  | PASS      |
| S4  | 1.00  | 0.510 | 0.7  | PASS      |
| S5  | 0.50  | 0.287 | 2.0  | PASS      |
| S6  | 0.25  | 0.132 | N/A  | PASS      |
| S7  | 0.125 | 0.096 | N/A  | PASS      |
| S8  | 0.00  | 0.051 | N/A  | PASS      |

Curve Details  
 Intercept: 0.7725 Slope: 0.7683  
 Correlation: 0.93 PASS

View test Quality Control.

| Subject ID | NI   | TB Ag | Mitogen | TB Ag/NI | Mitogen/NI | Result        |
|------------|------|-------|---------|----------|------------|---------------|
| A0090001   | 0.08 | 0.14  | > 10    | 0.06     | > 10       | NEGATIVE      |
| A0090002   | 7.23 | 0.19  | > 10    | 0.96     | 4.25       | NEGATIVE      |
| A0090003   | 0.06 | 0.31  | 0.87    | 0.25     | 0.01       | NEGATIVE      |
| A0090004   | 0.06 | 0.32  | 4.17    | 0.36     | 4.11       | NEGATIVE      |
| A0090005   | 0.06 | 0.10  | > 10    | 0.04     | > 10       | NEGATIVE      |
| A0090006   | 0.08 | 0.27  | > 10    | 0.19     | > 10       | NEGATIVE      |
| A0090007   | 0.09 | 0.13  | 0.39    | 0.04     | 0.29       | INDETERMINATE |
| A0090008   | 2.06 | 2.62  | 2.24    | 0.56     | 0.18       | POSITIVE      |
| A0090013   | 0.09 | 0.46  | 0.26    | 0.37     | 0.17       | POSITIVE      |
| A0090014   | 0.14 | 0.83  | 0.46    | 0.69     | 0.32       | POSITIVE      |

View test results.

Print test results, or save as PDF file.

Save File  
 Results Export  
 Data Export  
 Load File

Further options.

For comprehensive instructions, please refer to the Analysis Software Instructional Guide

North America / South America  
 Email: customer.service@cellestis.com  
 Fax: +1 661 775 7479  
 Tel: +1 661 775 7480 (outside USA)  
 Toll free: 800 519 4627 (USA only)

Europe / Middle East / Africa  
 Email: europe@cellestis.com  
 Fax: +49 6151 428 59 110  
 Tel: +49 6151 428 59 0

Asia / Oceania  
 Email: quantiferon@cellestis.com  
 Fax: +61 3 9571 3544  
 Tel: +61 3 9571 3500