



Cellestis International Pty Ltd
1046A Dandenong Rd
Carnegie Victoria
Australia 3163
Ph: 61 3 9571 3500
Fax: 61 3 9571 3544
quantiferon@cellestis.com
www.cellestis.com
ACN 100 319 122

CERTIFICATE OF ANALYSIS

QuantiFERON® - TB Gold: Mitogen tubes (Cat No.: T0593-0201)
Lot Number: 0593-60171
Date of Expiry: 2010-10

| | |
|-----------------------------|---|
| DESCRIPTION: | The Mitogen blood collection tubes are for use with the QuantiFERON®-TB Gold assay, which is a test of Cell Mediated Immunity for the detection of tuberculosis (TB) infection in humans. This simple and convenient blood test is designed to identify individuals who are infected with <i>Mycobacterium tuberculosis</i> complex (MTB) organisms. QuantiFERON®-TB Gold also provides an indication of the patient's general immune competence. |
| COMPOSITION: | The Mitogen tubes contain Phytohaemagglutinin-P (PHA-P) and heparin, which have been dried to the inner wall of the tubes. |
| PRESERVATIVE: | Not applicable |
| PRECAUTIONS: | Not applicable |
| STORAGE TEMPERATURE: | 4°C to 25°C |



Cellestis International Pty Ltd
1046A Dandenong Rd
Carnegie Victoria
Australia 3163
Ph: 61 3 9571 3500
Fax: 61 3 9571 3544
quantiferon@cellestis.com
www.cellestis.com
ACN 100 319 122

CERTIFICATE OF ANALYSIS

QuantiFERON® - TB Gold: Mitogen tubes (Cat No.: T0593-0201)

Lot Number: 0593-60171

Date of Expiry: 2010-10

QUALITY CONTROL TESTING:

Mitogen tubes

| Test | Criteria | Result |
|--------------------|--|--------|
| Visual examination | Clear coating visible on inside of tubes | Clear |
| Identity Test | Pass / Fail | Pass |
| Endotoxin Test | Pass / Fail | Pass |

NUMBER OF COMPONENTS: Components and Instruction leaflet

| Component | Quantity per carton | Lot / Version Number | Expiry Date |
|---------------------|---------------------|----------------------|----------------|
| Mitogen tubes | 100 | A0907079 | 2010-10 |
| Instruction Leaflet | 1 | L05900201K | Not applicable |

This certifies that the results obtained for this batch of QuantiFERON®-TB Gold tubes conform to the quality specifications set by Cellestis Quality Control.

Dr. Stephen Jones
Chief Technical Officer, Cellestis Limited

Date: 20 August 2009